





1000 ng/mL rabbit PV25 in:

- Potassium-phosphate pH 6.5 0.1 mol/L
- Bovine serum albumin (BSA) 5.0 g/L
- Kathon MW/WT 0.2 ml/L

*The PV25 antibody is ONLY added to the solution directly before use at a concentration of 1ug/ml!*

*Stability: at least 6 months at 2-8°C*

#### **REAGENT 4. Test dilution buffer**

- Tris/HCl pH 7.5 0.2 mol/L
- Bovine serum albumin (BSA) 5.0 g/L
- Kathon MW/WT 0.2 ml/L
- Phenol 1.0 g/L

*Preparation: add 6.05 g of Tris/HCl, 0.1 ml kathon MW/WT, 0.5 g phenol and 2.5 BSA in 400 ml ddH<sub>2</sub>O. Adjust pH to 7.5 with 4 N HCl, and complete the volume to 500 ml with ddH<sub>2</sub>O.*

*Stability: at least 6 months at 2-8°C*

#### **REAGENT 5. Tetramethylbenzidin-H<sub>2</sub>O<sub>2</sub> solution (TMB)**

- 3,3',5,5'-tetramethylbenzidine 20 mmol/L
- Aceton 100 ml/L
- Ethanol 900 ml/L
- H<sub>2</sub>O<sub>2</sub> 50 mmol/L

*Preparation: dissolve 240 mg TMB in 5 ml aceton. Further add 45 ml ethanol and then 300 µl of 30% H<sub>2</sub>O<sub>2</sub> solution.*

*Stability: at least 6 month at 15-25°, in a well closed, light-protected glass container*

#### **REAGENT 6. Substrate buffer**

- Potassium-citrate buffer pH 4.1 30 mmol/L
- Kathon MW/WT 0.2 ml/L

*Preparation: dissolve 6.3 g citric acid monohydrate in about 800 ml ddH<sub>2</sub>O.*

*Adjust pH to 4.1 with 4 N KOH. Complete the volume to 1000 ml with ddH<sub>2</sub>O and add 200 µl kathon MW/WT*

*Stability: Few months at 15-25°C*

### **Other solutions**

#### **Sensitization buffer (SOLUTION #3.1.)**

- NaHCO<sub>3</sub> pH 9.5 100 mmol/L

*Preparation: dissolve 4.2 g NaHCO<sub>3</sub> in 400 ml ddH<sub>2</sub>O. Adjust pH to 9.5 with 4 N NaOH. Complete the volume to 500 ml with ddH<sub>2</sub>O*

*Stability: at least 6 months at 2-8°C*



### **Blocking buffer (SOLUTION #3.2.)**

- Tris/HCl pH 7.5 200 mmol/L
- Bovine serum albumin (BSA) 10 g/L
- Kathon MW/WT 0.2 ml/L

*Preparation: dissolve 12.1 g of Tris-HCl and 0.1 ml Kathon MW/WT in 400 ml ddH<sub>2</sub>O. Adjust pH to 7.5 with 4 N HCl. Further add 5 g BSA. Complete the volume to 500 ml with ddH<sub>2</sub>O*

*Stability: at least 6 months at 2-8°C*

### **Stop buffer (SOLUTION #3.3.)**

- Sulfuric acid 1.0 mol/L

*Preparation: add 28 ml sulfuric acid to 500 ml ddH<sub>2</sub>O*

→ Warning: do not add water in acid

### **Washing buffer (SOLUTION #3.4.)**

- Tween-20 0.5 ml/L

*Preparation: dilute 500 µl Tween-20 in 1000 ml ddH<sub>2</sub>O*

### **Other materials**

**Micro-well plates:** polystyrol, flat bottom, 96-wells, NUNC-immunoplate, MAXISORP, 4-42 404

**Multichannel photometer:** filter 450 nm

### **ELISA running**

#### ***Sensitization of micro-well plate with rabbit PV25 antibody***

20 ml sensitization solution is required per micro-well plate

- A corresponding amount of antibody solution (**REAGENT 1**) is diluted 1:1000 in sensitization buffer (**SOLUTION #3.1.**): for example 20 µl **REAGENT 1** in 20 ml **SOLUTION #3.1.**  
→ rabbit PV25 antibody (end concentration 1:5'000 - 1:20'000) in 100 mmol/L NaHCO<sub>3</sub> pH 9.5
- Pipette 200 µl of this solution per well
- Incubate 16-24 hrs at 15-25°C
- Remove sensitization solution and wash 3x with ddH<sub>2</sub>O
- Pipette 200 µl of blocking buffer (**SOLUTION #3.2.**) per well
- Close the plate with parafilm and incubate at least one day at 15-25°C

The micro-well plates sensitized with PV25 antibody can be stored (under humidified conditions at 2-8°C) in the blocking buffer until 2 months without significant loss of immunological binding capacity

#### ***Preparation of solutions for standard curve and loading of samples***

- Dilute Parvalbumin standard solution 1000 ng/ml (**REAGENT 3**) in test dilution buffer (**REAGENT 4**)  
→ 1:200 for 5.0 ng/ml  
→ 1:400 for 2.5 ng/L



- 1:800 for 1.25 ng/L
  - The null-value is given by the test dilution buffer alone
  - h) Prepare mouse PV235 antibody test solution: Dilute mouse PV235 antibody solution (**REAGENT 2**) 1:500 in test dilution buffer (**REAGENT 4**) (mouse PV235 antibody end concentration 1:5'000 - 1:10'000)
- 5 ml of this preparation is required per micro-well plate (5 ml **REAGENT 4** and 5 µl **REAGENT 2**)
- i) Remove blocking buffer from the wells
  - j) Pipette 200 µl of standard curve solutions and samples (user-defined dilution with **REAGENT 4**) in duplicates in each well, and add 50 µl of PV235 antibody test solution prepared in step h)
  - k) Close the plate and incubate it for 16-24 hrs at 2-8°C (possible to do it over the week end)
  - l) Empty the wells and wash 3x with washing buffer (**SOLUTION #3.4.**)
  - m) Incubation with 200 µl of rabbit anti-mouse-HRP (DAKO P0260, 1:1000 in **REAGENT 4**) for 4-5 hrs at room temperature.
  - n) Empty the wells and wash 3x with washing buffer (**SOLUTION #3.4.**)
  - o) Wash 2-3x with dd. H<sub>2</sub>O.

**Enzyme reaction**

- p) Mix 1 volume of TMB-H<sub>2</sub>O<sub>2</sub> solution (**REAGENT 5**) with 20 volumes of substrate buffer (**REAGENT 6**)
- 21 ml of this solution is required per micro-well plate
- This solution should not be vortexed (could cause precipitation) and has to be prepared just before use (stability 1 hr at 15-25°C)
- q) Pipette 200 µl of this TMB-H<sub>2</sub>O<sub>2</sub> solution per well
  - r) Incubate 10 minutes at 15-25°C
  - s) Stop the reaction by adding 100 µl stop buffer (**SOLUTION #3.3.**)
  - t) Within 1 hr, analyze the plate with the photometer